Conservar entre: +2+8°C.

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Cod. QC008 CK-MB Control 1 x 3 mL.

CK-MB CONTROL SERUM

IVD

PRODUCT CHARACTERISTICS

The CK-MB control is a human lyophilised serum. Is intended for control of accuracy for use with manual and automated analytical procedures.

REAGENTS

Human serum and human derivated isoenzyme, with biological additives and bacteriostatic agents.

The concentration is lot-specific and its level is in the pathological range.

PRECAUTIONS

Components from human origin have been tested and found to be negative for the presence of HBsAg and HCV, and of antibody to HIV (1/2). However handle cautiously as potentially infectious.

PREPARATION

1. Remove vials from the refrigerator and allow to warm to room temperature for 15-20 minutes.
2. Open the bottle carefully, to avoid loss of lyophilized.
3. Reconstitute the lyophilised vial with 3 ml. of distilled water.
4. Mix thoroughly, avoiding foam forming. Wait for 20 minutes before use.

STORAGE AND STABILITY

The control serum is stable until the expiration date on the label when stored tightly closed at 2-8° C and contaminations are prevented during their use.

The stability after reconstitution is:

- 5 days at +2 to +8° C.

Do not use reagents over the expiration date.

Do not use the product if there is visible evidence of biological, chemical or physical deterioration.

Signs of Reagent deterioration:

- Presence of particles and turbidity.

TRACEABILITY¹

Concentration value for CK-MB CTRL is traceable to IFCC method.

EXPECTED VALUES

The values are indicated on the vial label.

LIMITATIONS

The control has been validated with CHEMELEX CK-MB LQ reagent. Users should verify its suitability with other reagents.

BIBLIOGRAPHY

1. Schumann, G. et al., IFCC Primary Reference Procedures for the Measurement of Catalytic Activity Concentrations of Enzymes at 37°C, Clin. Chem. Lab. Med. (2002), 40, 635.

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