Store at / Conservar a: +2+8°C.

*Presentación / Presentación:*  
Cod. TL012 CONT: 4 x 1 mL.Procedure / Procedimiento**Quantitative determination of ASO, CRP, RF.****Only for in vitro use in clinical laboratory (IVD)****CHARACTERISTICS OF PRODUCT**

The control is a human serum used to evaluate the precision and accuracy of ASO, CRP and RF determinations in human serum by turbidimetry method.

**REAGENTS COMPOSITION**

ASO/CRP/RF Human serum. With a low concentration of ASO, CONTROL L CRP and RF.

**PRECAUTIONS**

Components from human origin have been tested and found negative for the presence of HBsAg, HCV and antibody to HIV (1/2). However handle cautiously as potentially infectious.

**REAGENT PREPARATION AND STABILITY**

Reconstitute (→) accurately with 1.0 mL of distilled water. Mix thoroughly, avoiding foam forming. Bring at room temperature for about 10 minutes before use.

The control is stable until the expiration date on the label when stored tightly closed at 2-8°C and contaminations are prevented during their use. Do not use reagents over the expiration date.

After reconstitution, stable for 10 days at 2-8°C or 6 months at -20°C.

**CALIBRATION**

The concentration values of the Control have been standardized against the Reference Material ERM-DA 474/IFCC for CRP, International Reference Standard from NIBSC 64/002 for RF and International Standard from NIBSC ASO for ASO.

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